

See It. Believe It.

HydroMARK[®]

Breast Biopsy Site Marker

- HydroMARK is more than 90% water after hydration.
- Histology showed no significant migration.
- HydroMARK is a polyethylene glycol-based hydrogel, which contains a central marker of titanium or stainless steel.
- The polymer does not contain any collagen. It is biocompatible and completely resorbable.
- Marker expands in length and diameter after it is deployed, anchoring it within the biopsy cavity.
- HydroMARK has a 3 year shelf life.

HydroMARK Polymer Transition




Before




After

Pictures not to scale.

MRI



MRI




At Deployment



Ultrasound



Mammography



"We are extremely impressed with the HydroMARK markers. Not only are they well seen in all modalities, they are easy to use and priced competitively. We do virtually all our wire localizations under ultrasound now because the HydroMARK marker is so visible. We have a very active breast MRI practice and the clip jumps out at you!"

*Robin Shermis, MD
Toledo Hospital Breast Care Center*

Biopsy[®]
Sciences

4900 Creekside Drive, Suite C • Clearwater, FL 33760 USA
www.HydroMARK.com • www.BiopsySciences.com



HydroMARK®

Breast Biopsy Site Marker

See It. Believe It.

PRODUCT LIST

Biopsy®
Sciences

MARKER SHAPES

T1	T3	T4
Titanium or Stainless Steel	Titanium or Stainless Steel	Titanium

ULTRASOUND GUIDED BIOPSY

Compatibility	Delivery System	Coil	Reference #	Qty/Box	Shape
Ultrasound Guided Biopsy 15 GA	Rigid/Sharp Delivery System	Titanium	4010-02-15-T1	10	1
			4010-02-15-T3	10	3
			4010-02-15-T4	10	4
		SST	4010-02-15-S1	10	1
			4010-02-15-S3	10	3
Ultrasound Guided Biopsy 18 GA	Rigid/Sharp Delivery System	Titanium	4010-02-18-T3	10	3

MAMMOTOME® PROBE

Compatibility	Delivery System	Coil	Reference #	Qty/Box	Shape
Mammotome® 11 GA Biopsy Probe	Flexible Delivery System	Titanium	4010-01-11-T1	10	1
			4010-01-11-T3	10	3
			4010-01-11-T4	10	4
Mammotome® 8 GA Biopsy Probe & 8 GA MR Mammotome® Targeting Set	Flexible Delivery System	Titanium	4010-01-08-T1	10	1
			4010-01-08-T3	10	3
			4010-01-08-T4	10	4
		SST	4010-01-08-S1	10	1


Eviva® HANDPIECE*

Compatibility	Delivery System	Coil	Reference #	Qty/Box	Shape
Eviva® 9 GA	Rigid Delivery System End Deploy	Titanium	4010-04-09-T1	10	1
			4010-04-09-T3	10	3
			4010-04-09-T4	10	4
		SST	4010-04-09-S3	10	3

* Each Eviva device is packaged with adapters so that it may be used with either standard or petite probes.

ATEC® HANDPIECE INTRODUCERS

Compatibility	Delivery System	Coil	Reference #	Qty/Box	Shape
ATEC® 9 GA Biopsy Handpiece (Outer Cannula Method Only)	Flexible Delivery System	Titanium	4010-01-11-T1	10	1
			4010-01-11-T3	10	3
			4010-01-11-T4	10	4
ATEC® 12 GA Biopsy Handpiece (Outer Cannula Method Only)	Short Rigid Delivery System	Titanium	4010-03-15-T1-SHORT	10	1
			4010-03-15-T3-SHORT	10	3




ATEC® HANDPIECE INTRODUCERS

Compatibility	Delivery System	Coil	Reference #	Qty/Box	Shape
ATEC® 9 GA Biopsy Handpiece Introducer	Rigid Delivery System	Titanium	4010-03-09-T1	10	1
			4010-03-09-T3	10	3
ATEC® 12 GA Biopsy Handpiece Introducer	Rigid Delivery System	Titanium	4010-03-15-T1	10	1
			4010-03-15-T3	10	3

ATEC® MRI Biopsy Handpiece (Introducer Sheath Method Only)

Compatibility	Delivery System	Coil	Reference #	Qty/Box	Shape
ATEC® MRI Biopsy Handpiece (Introducer Sheath Method Only)	Flexible Delivery System	Titanium	4010-01-08-T1	10	1
			4010-01-08-T3	10	3
			4010-01-08-T4	10	4
		SST	4010-01-08-S1	10	1



EnCor® DIRECTIONAL VACUUM-ASSISTED BIOPSY DEVICES

Compatibility	Delivery System	Coil	Reference #	Qty/Box	Shape
EnCor® 10 GA	Rigid Delivery System	Titanium	4010-03-09-T1	10	1
			4010-03-09-T3	10	3
EnCor® 12 GA	Rigid Delivery System	Titanium	4010-03-15-T1	10	1
			4010-03-15-T3	10	3

EnCor® DIRECTIONAL VACUUM-ASSISTED BIOPSY DEVICES

Compatibility	Delivery System	Coil	Reference #	Qty/Box	Shape
EnCor® MRI Biopsy Device	Flexible Delivery System	Titanium	4010-01-08-T1	10	1
			4010-01-08-T3	10	3
			4010-01-08-T4	10	4
		SST	4010-01-08-S1	10	1

Vacora® VACUUM-ASSISTED BIOPSY DEVICE

Compatibility	Delivery System	Coil	Reference #	Qty/Box	Shape
Vacora® Biopsy Probe 118mm MR Coaxial Cannula	Flexible Delivery System	Titanium	4010-01-08-T1	10	1
			4010-01-08-T3	10	3
			4010-01-08-T4	10	4
		SST	4010-01-08-S1	10	1

BIOPSY SCIENCES®

4900 Creekside Drive, Suite C • Clearwater, FL 33760 USA • Tel: 727-290-9825 • Fax: 727-290-9828
www.HydroMARK.com • www.BiopsySciences.com

Patent numbers 5410016, 6060582, 6083524, 6162241, 6177095, 6605294, 6790185 and patents pending. Mammotome® is a registered trademark of Devicor Medical Products, Inc. Eviva® and ATEC® are registered trademarks of Hologic, Inc. Encor® and Vacora® are registered trademarks of C.R. Bard, Inc. or affiliates in the United States and/or other countries.

