

Mammotome Markers




MammoMARK/CorMARK

Breast Biopsy Site Identifier




Three Distinct Shapes

BOWTIE



TRIPLE TWIST




U-SHAPE




Available in Titanium (nickel free) with visibility in Ultrasound, X-ray and MRI.

Highly Visible Under Ultrasound

MammoMARK/CorMARK is highly echogenic which allows for precise localization under ultrasound.¹




Immediately after deployment



1. Corsi F, Sorrentino L, Bossi D, Sartani A, Foschi D. Preoperative localization and surgical margins in conservative breast surgery. *Int J Surg Oncol*. 2013;2013:793819. doi:10.1155/2013/793819

Anti-Migration Technology


The collagen marker plug increases in size by 300% within 60 seconds which enables precise placement in the tissue cavity.^{1,2}



Post mammography
60 seconds after deployment

Hemostatic Properties

The MammoMARK/CorMARK marker plug can assist in fluid management by quickly absorbing surrounding fluid.²



Pre-deployment

Post deployment

Compatibility

Core Needle in Ultrasound Procedures

Direct Puncture - 14G

Vacuum-Assisted Breast Biopsy Probe*

Mammotome	Mammotome revolve - 8G & 10G
Hologic	Eviva [®] , Brevera [®] - 9G & 12G ATEC [®] - 9G & MRI
Bard (BD)	EnCor Enspire™, EnCor Ultra™ - 7G & 10G

*See the Mammotome Marker compatibility brochure for more information

Please note: The CorMARK applicator is designed for use with core needles in ultrasound biopsies. The MammoMARK applicators are designed for use through the back of a vacuum-assisted biopsy probe.

Mammotome

www.mammotome.com
1-877-9-A-MAMMO

© 2020 Devicor Medical Products, Inc. All rights reserved. Devicor, Mammotome, Mammotome revolve, CorMARK and MammoMARK are registered trademarks of Devicor Medical Products, Inc., in the USA and optionally in other countries. Hologic, Eviva[®], Brevera[®] and ATEC[®] are registered trademarks of Hologic, Inc. Bard, EnCor Enspire™ and EnCor Ultra™ are the product names, trade names and/or trademarks of C.R. Bard, Inc., and its affiliates including BD. Other logos, product and/or company names might be trademarks of their respective owners. MDM# 200566 Rev 11/20