



Monitores

Monitores | NDS

	Luz LED ultrabrillante	Luz LED		
	RADIANCE® ULTRA 27"	RADIANCE® G2 24"	RADIANCE® G2 HB 26"	RADIANCE® HD 19"
Señal de entrada	HD-SDI x 2, DVI-D, DVI-I, RGBS, YPbPr, S-Video, Composite, VGA	HD-SDI x 2, DVI-D, DVI-I, RGBS, YPbPr, S-Video, Composite, VGA	HD-SDI x 2, DVI-D, DVI-I, RGBS, YPbPr, S-Video, Composite, VGA	HD-SDI x 2, DVI-D, DVI-I, RGBS, YPbPr, S-Video, Composite, VGA
Señal de salida	HD-SDI, DVI, RGBS, YPbPr/VGA, S-Video, Composite	HD-SDI, DVI, RGBS, YPbPr/VGA, S-Video, Composite	HD-SDI, DVI, RGBS, YPbPr/VGA, S-Video, Composite	HD-SDI, DVI, RGBS, YPbPr/VGA, S-Video, Composite
Dimensiones (W x H x D)	(678 x 445 x 84 mm)	(597 x 401 x 100 mm)	(673 x 418 x 88 mm)	(465 x 400 x 98 mm)
Peso	8.9 kg	7.1 kg	8.2 kg	6.8 kg