

maccura

SCIENCE & TECHNOLOGY
Serve Human Health



Maccura E2 Series Blood Glucose Test System



Traceability Ensures the Results Accurate

The results could be traced back to the reference methods recommended by JCTLM



Easy and Convenient to Operate

Accurate results in only 5 seconds
Automatic unloading of test strips
No coding
Sample volume: 0.5 μ L



Meet the Latest International Standard

Accuracy requirements of ISO 15197: 2013



Able to all Kinds of Complex Environment

Test temperature: 5 $^{\circ}$ C-45 $^{\circ}$ C
Humidity: 10%-90%
Altitude: 0-3463 meters

SPECIFICATION

Model	E 201	E 202	E 203	E 205
Phonetic Function	/	/	Yes	/
Data Transmission	/	/	/	Bluetooth
Coding Type	Manual coding	No coding	No coding	No coding
Test Method	Glucose Oxidase Method	Glucose Oxidase Method	Glucose Oxidase Method	Glucose Oxidase Method
Sample Type	Capillary whole blood	Capillary whole blood	Capillary whole blood	Capillary whole blood
Sample Volume	0.5μL	0.5μL	0.5μL	0.5μL
Test Time	5 seconds	5 seconds	5 seconds	5 seconds
Test Range	1.1-33.3mmol/L	1.1-33.3mmol/L	1.1-33.3mmol/L	1.1-33.3mmol/L
HCT	25%-60%	25%-60%	25%-60%	25%-60%
Temperature	5°C-45°C	5°C-45°C	5°C-45°C	5°C-45°C
Altitude	≤3463 meters	≤3463 meters	≤3463 meters	≤3463 meters
Humidity	10%-90%	10%-90%	10%-90%	10%-90%
Memory Capacity	300	300	300	300
Strips Unloading	Automatic	Automatic	Automatic	Automatic

MACCURA BIOTECHNOLOGY CO., LTD.

Add: 8#, 2nd Anhe Road, Hi-tech Zone, 611731 Chengdu
Tel: +86 28 8173 1888
E-mail: maccura@maccura.com www.maccura.com



WEBSITE

WECHAT