



# URINALYSIS

Efficient Intelligent Sensitive

SCIENCE & TECHNOLOGY  
*Serve Human Health*

**maccura**

# Efficient Intelligent Sensitive

Efficient test speed , the throughput can up to 600 T/H  
Automatically correct the effect of uric acid alkalinity on specific gravity test results  
Built-in temperature detection system can correct the results according to the temperature  
New type of light sensor, high sensitivity, strong anti-interference


## Parameters

Test Items: (URO)/(BIL)/(KET)/(BLD)/(PRO)/(NIT)/(LEU/GLU)/(pH)/(SG)/(VC)/(MA)/(Cr)/(Ca)/(SG)  
Compatible Reagents : 8-14 options available  
Principles: reverse color ratio and color analysis  
Input mode:U 390 Touch Input (LCD)  
Data Storage Capacity: 2000 test reports  
Report unit: support international system of units, conventional system of units, symbol system  
Sensitivity: Adjust the software Settings, after each project was opened  
Light source: high brightness cold light source, long service life

## U 390 Semi-Automatic Urine Analyzer



## U 350 Semi-Automatic Urine Analyzer



Efficient test speed , the throughput can up to 120 T/H  
Automatically correct the effect of uric acid alkalinity on specific gravity test results  
New type of light sensor, high sensitivity, strong anti-interference

## Parameters

Test Items: (URO)/(BIL)/(KET)/(BLD)/(PRO)/(NIT)/(LEU)/(GLU)/(pH)/(SG)/(VC)/(MA)/(Cr)/(Ca)/(SG)  
Compatible Reagents: 8-14 options available  
Principles: reverse color ratio and color analysis  
Throughput: 120T/H (fast mode), 60T/H (slow mode), manual trigger  
Data Storage Capacity : 10,000 test reports  
Report unit: support international system of units, conventional system of units, symbol system  
Sensitivity: Adjust the software Settings, after each project was opened  
Light source: high brightness cold light source, long service life

**NOTICE:** This brochure is for reference purpose only and subject to change without notice.

MACCURA BIOTECHNOLOGY CO., LTD.

Add: 8#, 2nd Anhe Road, Hi-tech Zone, 611731 Chengdu  
Tel: +86 28 8173 1888  
E-mail: [maccura@maccura.com](mailto:maccura@maccura.com) [www.maccura.com](http://www.maccura.com)



WEBSITE



WECHAT

